

SURGICAL TECHNIQUE GUIDE

FREEFIX[®]

ulnar shortening system



 **skeletal dynamics**[®]
UNDERSTANDING THE UPPER EXTREMITY

As described by:
Jorge L. Orbay, M.D.
Miami Hand & Upper Extremity Institute

FREEFIX[®]

ulnar shortening system

Description

The Skeletal Dynamics Ulnar Shortening Plating System consists of the following:

1. Distal Ulnar Shortening Plates
2. Proximal Ulnar Shortening Plates
3. Screws, K-wires and Specialized Instrumentation

Both Distal and Proximal Ulnar Shortening Plates are titanium alloy plates and designed for fracture fixation, osteotomies, and non-unions of the ulna.

Distal Ulnar Shortening Plates

The Distal Ulnar Shortening Plates are available in two lengths, left and right configurations.

For surgical instructions for the Distal Ulna Shortening Plate see page 4.

Proximal Ulnar Shortening Plates

The Proximal Ulnar Shortening Plates are available in two lengths.


For surgical instructions for the Proximal Ulna Shortening Plate see page 11.

Indications for Use

The Skeletal Dynamics Ulnar Shortening Plating System is indicated for fractures and osteotomies, in particular for the ulna.


DISTAL

Right and left specific
Two lengths



PROXIMAL

Non side specific
Two lengths



1

EXPOSURE



Make an 8cm to 10cm longitudinal incision ulnar distally between the extensor carpi ulnaris (ECU) and flexor carpi ulnaris (FCU) compartment.

Warning:

Retract and protect the tendons and dorsal branch of the ulnar nerve.

2

PLATE POSITIONING



Determine proper distal plate placement using fluoroscopy. Mount the k-wire AIMing guide in the most distal pre-loaded drill guide using the torque limiting driver and handle. Estimate the desired plate position and insert the 1.6mm k-wire. Using fluoroscopy, confirm the k-wire courses just proximal to the distal ulnar articular surface. If the plate is misaligned, remove the distal k-wire and reposition the plate, confirming with fluoroscopy. Bone

holding forceps have been provided to assist in fixation.

Note:

Utilize the PROTEAN® Bending Pliers to adjust the inclination of the neck of the plate as needed.

Warning:

Caution should be taken when contouring plates. Bending the plates may weaken or break the plates.



PDG-AIM-015: AIMing Guides, 1.5mm



PRT-BND-PLL:
PROTEAN® Plate
Bending Pliers, Large

KWIR-SD-16127: K-Wire, Single Diamond, 1.6mm x 127mm

DISTAL FREEFIX® SLOT

3

Place the chosen drill guide in the most proximal position of the distal FreeFix® slot and drill with the 2.7mm drill. Measure the screw length and utilize the T-10 driver to place the selected 3.5mm cortical non-locking screw.



Warning:

At this time, do not use a cortical locking screw in the FreeFix® slot to prevent deformation of the cutting guide rails. If desired, a cortical locking screw may replace the non-locking screw after the final osteotomy has been achieved.

If replacing a cortical non-locking screw with a locking screw in a later step, use the thread-in drill guide.



TPDG-FF-35: Thread-in Drill Guide, FreeFix®, 3.5mm



DRVR-UQC-T10: Driver, Universal Quick Connect, T10



DRLL-SSC-27040: Drill, 2.7mm x 40mm



DPGA-FF-050: FreeFix® Depth Gauge




TPDG-FF-USP: Drill Guide, FreeFix®, Ulnar Shortening

PROXIMAL FREEFIX® SLOT

4

Place the chosen drill guide in the most proximal position of the proximal FreeFix® slot and drill with the 2.7mm drill. Measure the screw length and utilize the T-10 driver to place the selected 3.5mm cortical non-locking screw. Ensure the screw is 2mm longer than measured to prevent the screw from stripping during compression (see step 12).



Warning:

At this time, do not use a cortical locking screw in the FreeFix® slot to prevent deformation of the cutting guide rails. If desired, a cortical locking screw may replace the non-locking screw after the final osteotomy has been achieved.

If replacing a cortical non-locking screw with a locking screw in a later step, use the thread-in drill guide.

5

SCREW PLACEMENT



Use the 2.0mm drill bit to drill through the three pre-loaded drill guides on the head of the plate. Measure each screw length using the depth gauge or measurement lines on the drill bit, taking note of the appropriate scale, and utilize the peg driver to insert a 2.7mm non-locking screw or 2.3mm locking screw of selected length in each hole. A thread-in drill guide is available if needed.

Warning:

Prevent excessive screw length as this will cause the screws to violate the distal radioulnar joint (DRUJ) space.



DRVR-AOS-S20: Driver, Peg, Torque Limiting




TPDG-THD-DG20: Thread-in Drill Guide, 2.0mm



DRLL-SSC-20040: Drill, Solid Side Cutting, 2.0mm x 40mm

6

DISTAL SCREW



Remove the k-wire from the head of the plate, remove AIMing guide, then drill using the 2.0mm drill bit. Measure the screw length using the depth gauge, taking note of the appropriate scale, and insert a 2.3mm locking screw or 2.7mm non-locking screw of selected length in the distal hole with the peg driver. Remove all guides.



Warning:

Prevent excessive screw length as this will cause the screws to violate the distal radioulnar joint (DRUJ) space.

Attach the appropriate distal cutting guide (based on desired osteotomy approach) to the plate using the rails along the sides of the plate. Push the cutting guide until it reaches the end of the rail. It stops just distal to the leading edge of the thick, distal line.

To fix the cutting guide, either hold it in place by pushing it against the end of the rail with a thumb or use two 1.6mm k-wires through the k-wire holes on the side of the cutting guide.

Insert the k-wires in two different columns to avoid k-wire convergence.



DISTAL CUT

8

Select the suggested saw blade by referencing the chart on page 23.

Use an oscillating saw and the selected saw blade to resect the ulna through the guiding slot on the distal cutting guide arm. Once the ulna has been resected all the way through, remove the saw blade and any k-wires from the bone.

Warning:

Protect soft tissues while cutting with the oscillating saw using the Hohmann Retractors provided in the set.



INST-HR-DBL: Instrument, Hohmann Retractor, Double Ended




INST-MHR-STD: Instrument, Mini-Hohmann Retractor, Standard (optional)



CUTG-DIST-R/L: Ulnar Shortening, Cutting Guide, Distal

9

ADJUST FOR LENGTH




Adjust the cutting guide proximally until the determined osteotomy length is achieved. Each laser-etched line on the plate corresponds to an additional osteotomy length of 2mm.

To secure the cutting guide, use two 1.6mm k-wires through the k-wire holes on the side of the cutting guide. Note the k-wire placement; see step 7 for details.

10

PROXIMAL CUT



Use an oscillating saw and the saw blade selected in step 8 to resect the ulna through the guiding slot on the distal cutting guide arm. Once the ulna has been resected all the way through, remove the saw and any k-wires from the bone.


Warning:

Protect soft tissues while cutting with the oscillating saw using the Hohmann Retractors provided in the set.

BONE REMOVAL


11

Fully remove the resected bone wafer from the osteotomy site. A Cottle Osteotome has been provided in the set to facilitate removal.




OSTEOTOMY REDUCTION

12



Loosen the cortical non-locking screws by approximately one half of a turn or until there is little friction for the plate to translate upon the bone in the distal and proximal FreeFix® slots. Place the appropriate tip of the Compression Forceps into the recess at the center of the plate and the other tip behind the head of the screw in the proximal slot. Using the Compression Forceps, reduce and compress the osteotomy site until adequate compression is achieved. Then tighten the distal and proximal cortical non-locking screws.



13 FINAL CORTICAL SCREW



Use the 2.7mm drill bit (and thread-in drill if selecting a cortical locking screw) to drill for the final 3.5mm cortical screw. Measure the appropriate length using the depth gauge and place the screw into the proximal hole of the plate using the T-10 driver.

Note:

Additional dynamic compression may be achieved across the osteotomy site by either drilling the final screw eccentrically in the hole, or placing a cortical screw (either locking or non-locking) on the proximal edge of either of the FreeFix® slots. Be sure to loosen the screws in the distal and proximal FreeFix® slots prior to achieving dynamic compression.

14 CLOSURE



Check the plate and screw trajectories under fluoroscopy for adequate placement and osteotomy compression. If satisfied, close the site in the usual fashion.

Make a 10cm to 12cm longitudinal incision along the subcutaneous border of the distal to mid-ulna. Determine proper plate placement using fluoroscopy. Plate should lie proximal to the distal band and distal to the central band of the interosseous ligament. Bone holding forceps have been provided to assist in provisional fixation.

The osteotomy should be placed on the ulna 3/4 of way distal on the forearm. This will preserve the integrity of the central band of the interosseous ligament. The plate can be placed anterior or posterior on the ulna per surgeon's preference.



Use the 2.7mm drill bit (and thread-in drill guide if selecting a cortical locking screw) to drill through the three distal holes on the plate. Measure each screw length using the appropriate scale of the depth gauge and insert a 3.5mm cortical non-locking or 3.5mm cortical locking screw of selected length in each hole using the T-10 driver.



DRVR-UQC-T10: Driver, Universal Quick Connect, T10



DRLL-SSC-27040: Drill, 2.7mm x 40mm



DPGA-UNV-030: Depth Gauge, Universal, 30mm



TPDG-FF-35: Thread-in Drill Guide, FreeFix®, 3.5mm



3

DISTAL FREEFIX® SLOT



Place the chosen drill guide in the most proximal position of the distal FreeFix® slot and drill with the 2.7mm drill. Measure the screw length and utilize the T-10 driver to place the selected 3.5mm cortical non-locking screw.

Warning:

At this time do not use a cortical locking screw in the FreeFix® slot to prevent deformation of the cutting guide rails. If desired a cortical locking screw may replace the non-locking screw after the final osteotomy has been achieved.



If replacing a cortical non-locking screw with a locking screw in a later step, use the thread-in drill guide.



TPDG-FF-USP: Drill Guide, FreeFix®, Ulnar Shortening



DPGA-FF-050: FreeFix® Depth Gauge

4

PROXIMAL FREEFIX® SLOT



Place the chosen drill guide in the most proximal position of the proximal FreeFix® slot and drill with the 2.7mm drill. Measure the screw length and utilize the T-10 driver to place the selected 3.5mm cortical non-locking screw. Ensure the screw is 2mm longer than measured to prevent the screw from stripping during compression (see step 10).

Warning:

At this time do not use a cortical locking screw in the FreeFix® slot to prevent deformation of the cutting guide rails. If desired a cortical locking screw may replace the non-locking screw after the final osteotomy has been achieved.



If replacing a cortical non-locking screw with a locking screw in a later step, use the thread-in drill guide.



Check the placement of the plate using fluoroscopy. If the plate is misaligned, readjust as necessary.

CUTTING GUIDE

Attach the appropriate proximal cutting guide to the plate using the rails along the side of the plate. Both transverse and oblique cut options are available in left and right configurations depending on the preferred osteotomy approach. Push the cutting guide until it reaches the end of the rail. It stops just distal to the leading edge of the thick, distal line.

To secure the cutting guide, either hold it in place by pushing it against the end of the rail with a thumb or use two 1.6mm k-wires through the k-wire holes on the side of the cutting guide.

Insert the k-wires in two different columns to avoid k-wire convergence.



Correct

Correct

Incorrect

KWIR-SD-16127: K-Wire, Single Diamond, 1.6mm x 127mm


DISTAL CUT

Select the suggested saw blade by referencing the chart on page 23.


Use an oscillating saw and the selected saw blade to cut the ulna through the guiding slot on the proximal cutting guide arm. Once the ulna has been cut all the way through, remove the saw and any k-wires from the bone.

Warning:

Protect soft tissues while cutting with the oscillating saw using the Hohmann Retractors provided in the set.



CUTG-OPRX-R/L: Ulnar Shortening, Cutting Guide, Oblique, Proximal



CUTG-PRX-R/L: Ulnar Shortening, Cutting Guide, Proximal

7

ADJUST FOR LENGTH



Adjust the cutting guide proximally until the determined osteotomy length is achieved. Each laser etched line on the plate corresponds to an additional osteotomy length of 2mm.

To secure the cutting guide, use two 1.6mm k-wires through the k-wire holes on the side of the cutting guide. Note the k-wire placement; see step 5 for details.

8

PROXIMAL CUT



Use an oscillating saw and the saw blade selected in step 6 to cut the ulna through the guiding slot on the proximal cutting guide arm. Once the ulna has been cut all the way through, remove the saw and any k-wires from the bone.

Warning:

Protect soft tissues while cutting with the oscillating saw using the Hohmann Retractors provided in the set.




INST-HR-DBL: Instrument, Hohmann Retractor, Double Ended



INST-MHR-STD: Instrument, Mini-Hohmann Retractor, Standard (optional)



Fully remove the cut bone wafer from the osteotomy site. A Cottle Osteotome has been provided in the set to facilitate removal.



OSTEOTOMY REDUCTION

10

Loosen the cortical non-locking screws by approximately one half of a turn or until there is little friction for the plate to translate upon the bone in the distal and proximal FreeFix® slots. Place the appropriate tip of the Compression Forceps into the recess at the center of the plate and the other tip behind the head of the screw in the proximal slot. Using the Compression Forceps, reduce and compress the osteotomy site until adequate compression is achieved. Then tighten the distal and proximal cortical non-locking screws.



11 OBLIQUE SCREW (OPTIONAL)



If an oblique osteotomy cut was chosen, an oblique lag screw may be placed across the osteotomy interface.

Use the 2.7mm drill bit and Oblique Compression Hole Guide to drill a pilot hole perpendicular to the plane of the osteotomy through both bone fragments. The Oblique Compression Hole Guide must be inserted at an angle and levered down to properly fit into the slot.

To ensure the screw will be properly lagged, overdrill the hole in the proximal osteotomy fragment only with the 3.5mm drill bit.




TPDG-CHG: Oblique Compression Hole Guide



DRLL-SSC-35040: Drill, 3.5mm x 40mm

12 OBLIQUE SCREW (OPTIONAL)



Measure the hole using the depth gauge and drive a 3.5mm cortical non-locking screw into the angled oblique hole of the plate. As the screw begins to achieve traction loosen the non-locking screws in the proximal and distal FreeFix® slots to compress along the osteotomy. Once the oblique screw is fully seated, tighten both screws in the proximal and distal FreeFix® slots.

To facilitate maintaining compression, utilize the compression forceps while driving the oblique screw into the bone.

FINAL CORTICAL SCREW 13

Use the 2.7mm drill bit (and thread-in drill guide if selecting a cortical locking screw) to drill through the most proximal hole for the final 3.5mm cortical non-locking screw. Measure the appropriate length using the depth gauge and insert the screw using the T-10 driver.

Note:

Additional dynamic compression may be achieved across the osteotomy site by either drilling the final screw eccentrically in the hole, or placing a cortical screw (either locking or non-locking) on the proximal edge of either of the FreeFix® slots. Be sure to loosen the screws in the distal and proximal FreeFix® slots prior to achieving dynamic compression.



CLOSURE 14

Check the plate and screw trajectories under fluoroscopy for adequate placement and osteotomy compression. If satisfied, close the site in the preferred manner.



1

SECOND HOLE IN DISTAL FREEFIX® SLOT



If the ulna requires more than 1cm of resection, proceed with the technique up to step 12, Osteotomy Reduction for the Distal Plate and to step 10, Osteotomy Reduction for the Proximal Plate.

With the ulna fully compressed after the first osteotomy, use the 2.7mm drill bit with thread-in drill guide to drill an additional hole in the most proximal position of the distal FreeFix® slot. Measure to confirm the appropriate screw length. If the measured length is the same

as the previously placed screw, the same screw may be used. Place the screw utilized during the initial osteotomy.

2

SECOND HOLE IN PROXIMAL FREEFIX® SLOT




Use the 2.7mm drill bit to drill an additional hole in the most proximal position of the proximal FreeFix® slot. If different screw lengths are needed, remove the screws initially placed in the slots and discard. Ensure the screw is 2mm longer than measured to prevent the screw from stripping during compression.

POSITIONING CUTTING GUIDE

3

Attach the distal cutting guide to the plate using the rails along the sides of the plate. Adjust the cutting guide proximally until the determined osteotomy length is achieved. Each laser-etched line on the plate corresponds to an additional osteotomy length of 2mm.

To secure the cutting guide, use two 1.6mm k-wires through the k-wire holes on the side of the cutting guide.



FINAL CUT


4

Use an oscillating saw to resect the ulna through the guiding slot on the distal cutting guide arm. Once the ulna has been resected all the way through, remove the saw and any k-wires from the bone.

Warning:

Protect soft tissues while cutting with the oscillating saw using the Hohmann Retractors provided in the set.

Fully remove the resected bone wafer from the osteotomy site. A Cottle Osteotome has been provided in the set to facilitate removal.



5

SECOND OSTEOTOMY REDUCTION




Loosen the cortical non-locking screws by approximately one half of a turn or until there is little friction for the plate to translate upon the bone in the distal and proximal FreeFix® slots.

Place the appropriate tip of the Compression Forceps into the recess at the center of the plate and the other tip behind the head of the screw in the proximal slot. Using the Compression Forceps, reduce and compress the osteotomy site until adequate compression is achieved. Then tighten the

distal and proximal cortical non-locking screws.

Continue the technique at step 13, Final Cortical Screw for the Distal Plate and step 11, Oblique Screw for the Proximal Plate.


INSTRUMENT TRAY (Standard Configuration)



Loc #	Catalog #	Description	Loc #	Catalog #	Description
1	USP-DRT USP-DLT	Ulnar Shortening Plate, Distal, Right Ulnar Shortening Plate, Distal, Left	11	TPDG-FF-USP	Drill Guide, FreeFix®, Ulnar Shortening
2	USP-DRS USP-DLS	Ulnar Shortening Plate, Distal, Short, Right Ulnar Shortening Plate, Distal, Short, Left	12	DPGA-FF-050	FreeFix® Depth Gauge
3	USP-PRX	Ulnar Shortening Plate, Proximal	13	DPGA-UNV-030	Depth Gauge, Universal, 30mm
4	USP-PRXL	Ulnar Shortening Plate, Proximal, Long	14	C-OST-20	Cottle Osteotome, 2.0mm
5	TPDG-THD-DG20 TPDG-FF-35	Thread-in Drill Guide, 2.0mm Thread-in Drill Guide, FreeFix®, 3.5mm	15	INST-MHR-STD INST-HR-DBL	Instrument, Mini-Hohmann Retractor, Standard (optional) Instrument, Hohmann Retractor, Double Ended
6	DRLL-SSC-20040 DRLL-SSC-27040 DRLL-SSC-35040	Drill, Solid Side Cutting, 2.0mm x 40mm Drill, 2.7mm x 40mm Drill, 3.5mm x 40mm	16	FRCP-BHM-RTC	Forceps, Bone Holding Medium, Ratcheting
7	HNDL-UQC-FXD	Handle, Universal Quick Connect, Fixed	17	FRCP-BHM-COMR FRCP-BHM-COML	Forceps, Bone Holding Medium, Compression, Right Forceps, Bone Holding Medium, Compression, Left
8	HNDL-SQC-FXD	Handle, Small QC, Fixed		PRT-BND-PLL	PROTEAN Plate Bending Pliers, Large *
9	DRVR-AOS-S20 DRVR-S20 DRVR-UQC-T10	Driver, Peg, Torque Limiting Driver, Peg Driver, Universal Quick Connect, T10			
10	TPDG-CHG	Oblique Compression Hole Guide			

* Not Pictured, stored in area below this layer.

INSTRUMENT TRAY (Standard Configuration)



Loc #	Catalog #	Description	Loc #	Catalog #	Description
18	CUTG-DIST-R	Ulnar Shortening, Cutting Guide, Distal, Right	26	COLS-35100-TS	Screw, Cortical Locking, 3.5mm x 10mm, Ti
19	CUTG-DIST-L	Ulnar Shortening, Cutting Guide, Distal, Left		COLS-35120-TS	Screw, Cortical Locking, 3.5mm x 12mm, Ti
20	CUTG-PRX-R	Ulnar Shortening, Cutting Guide, Proximal, Right		COLS-35140-TS	Screw, Cortical Locking, 3.5mm x 14mm, Ti
21	CUTG-PRX-L	Ulnar Shortening, Cutting Guide, Proximal, Left		COLS-35160-TS	Screw, Cortical Locking, 3.5mm x 16mm, Ti
22	CUTG-OPRX-R	Ulnar Shortening, Cutting Guide, Oblique, Proximal, Right		COLS-35180-TS	Screw, Cortical Locking, 3.5mm x 18mm, Ti
23	CUTG-OPRX-L	Ulnar Shortening, Cutting Guide, Oblique, Proximal, Left		COLS-35200-TS	Screw, Cortical Locking, 3.5mm x 20mm, Ti
24	TPNL-27100-TS	Threaded Peg, Non-Locking, 2.7mm x 10mm, Ti	27	PANL-35100-TS	Screw, Cortical Non Locking, 3.5mm x 10mm, Ti
	TPNL-27120-TS	Threaded Peg, Non-Locking, 2.7mm x 12mm, Ti		PANL-35120-TS	Screw, Cortical Non Locking, 3.5mm x 12mm, Ti
	TPNL-27140-TS	Threaded Peg, Non-Locking, 2.7mm x 14mm, Ti		PANL-35140-TS	Screw, Cortical Non Locking, 3.5mm x 14mm, Ti
	TPNL-27160-TS	Threaded Peg, Non-Locking, 2.7mm x 16mm, Ti		PANL-35160-TS	Screw, Cortical Non Locking, 3.5mm x 16mm, Ti
	TPNL-27180-TS	Threaded Peg, Non-Locking, 2.7mm x 18mm, Ti		PANL-35180-TS	Screw, Cortical Non Locking, 3.5mm x 18mm, Ti
	TPNL-27200-TS	Threaded Peg, Non-Locking, 2.7mm x 20mm, Ti		PANL-35200-TS	Screw, Cortical Non Locking, 3.5mm x 20mm, Ti
	TPNL-27220-TS	Threaded Peg, Non-Locking, 2.7mm x 22mm, Ti		PANL-35220-TS	Screw, Cortical Non Locking, 3.5mm x 22mm, Ti
25	TPLS-23100-TS	Threaded Peg, Locking, 2.3mm x 10mm, Ti	28	KWIR-SD-16127	K-Wire, Single Diamond, 1.6mm x 127mm
	TPLS-23120-TS	Threaded Peg, Locking, 2.3mm x 12mm, Ti	29	PDG-AIM-015	AIMing Guides, 1.5mm
	TPLS-23140-TS	Threaded Peg, Locking, 2.3mm x 14mm, Ti			
	TPLS-23160-TS	Threaded Peg, Locking, 2.3mm x 16mm, Ti			
	TPLS-23180-TS	Threaded Peg, Locking, 2.3mm x 18mm, Ti			
	TPLS-23200-TS	Threaded Peg, Locking, 2.3mm x 20mm, Ti			
	TPLS-23220-TS	Threaded Peg, Locking, 2.3mm x 22mm, Ti			

Ulnar Shortening Saw Blade Suggestions (not provided)

Stryker CORE/CORE 2/F1 System/SABO Saw Series	Conmed Linvatec MicroFree, MicroPower, MicroPower+	Conmed Linvatec Hall 50 Sagittal Saw	Synthes Small Battery System/TRS	Synthes Electric Pen Drive	Arthrex All Systems Small Bone Saw	SNN Dyonics Power Small	MicroAire Series 5000
2296-3-125 DISTAL	5023-249 (MicroAire ZS-125) DISTAL	5071-133 DISTAL	519.103S DISTAL	03.000.305S DISTAL	AR-300-402S DISTAL	72204242 (MicroAire ZS-125) DISTAL	ZS-125 DISTAL
31 x 9.8 x 0.40mm (MicroAire SP-125A)	31 x 9 x 0.40mm	25.5 x 9.5 x 0.40mm	25 x 10 x 0.60mm	31 x 10 x 0.38mm	25 x 9.4 x 0.6mm	31 x 9 x 0.40mm	31 x 9 x 0.40mm
2296-33-254 PROXIMAL	5023-133 (MicroAire ZS-333) PROXIMAL	5071-134 PROXIMAL	519.105S PROXIMAL	03.000.315S PROXIMAL	AR-300-400S PROXIMAL	72204243 (MicroAire ZS-333) PROXIMAL	ZS-333 PROXIMAL
41.6 x 19 x 0.51mm (MicroAire SP-410R)	41 x 14 x 0.40mm	41 x 9.0 x 0.40mm	50 x 20 x 0.60mm	43 x 20 x 0.60mm	41 x 14 x 0.6mm	41 x 14 x 0.40mm	41 x 14 x 0.40mm

